

DC-2000 CATHETERIZATION TABLE

Able supply to GE to work with GE-OEC 9600
C-arm since 1996GE-OEC 9800
C-armGE-OEC 9900 C-arm



Equip with various C-arm such as America GE-OEC,
SIEMENS, PHILIPS, GERNAMY ZIEHM and Chinese brand.

DC-2000 Catheterization table adopts the advanced electrical mechanical Transmission technology, the table top is imported carbon fibre, special for equip with kinds of mobile C-arm, small C-arm which both abroad and China. For the Micro-operation of intervention such as systemic vascular, ERCP, liver, tumor, angeitides, pain etc. This table is excellent technology, easy operate, elegant surface, advanced functions and superior quality. Had been supply for C-arm to GE, SIEMENS, PHILIPS, ZIEHM etc international brand since 1996, table is in good quality and welcome by user in China and abroad.

Supply tables for medical system of SIEMENS, GE, PHILIPS.

HOTBORN

The operating table adopts the imported micro-motor, with electric gear mesh mechanism structure instead of the traditional hydraulic transmission. Not need to change oil termly, and worry about the problem of leak oil. Product with fine steel can move equally and smoothly in any directions; also provides accurate positioning, trusty capability, and Durable lifetime.

The table top adopts the carbon fibre material, and 1100 mm length without metal. With low X-ray absorption, and gets high clear imaging picture. It extends the lifetime of Spheroid tube, battery and main parts of video machine, also reduces the radiation damage of X-ray to the patients and medical staff.



The big wheel especially use for medical treatment, the whole table is easily and convenient to be moved to clear the floor.

High precision imported beeline guide rail, ensure the smooth and balanced movement of operating table.

Base with the ABS industrial plastic, helps prevent glare, impact; and it is easily clear up and listerize.

Handle control is easily operating, and takes control of the horizontal move of table top; the brake system is safe and sensitive.

		Item	Data	
		Specification	Table top size(L×W×H)	
Length without metal			1100mm	
Length of imaging			1450mm	
Table top movement	Longitudinal		Head section end	Max.1700mm
			Foot section end	Max.810mm
			Travel	820mm
	Lateral		Travel	220mm
	Vertical		Min.distance from floor	795mm
			Travel	275mm
Noise			< 60dB	
Power		220V 50Hz 350VA		
Load		200kg		

STANDARD ACCESSORIES

Mattress	1 set
Arm support	1 set
Infusion pole	1 set
Handle controller	1 set

HOTBORN

Add: No.8 Keyuan Dong Wu Road , Nanning New& Hi-Tech Industrial Development Zone, Nanning , Guangxi , China
 Tel: +86-771-3947312 13877577356
 Fax: +86-771-3947312
 E-mail: NNHOTBO@VIP.163.COM
 Web: WWW.HOTBORN.COM

Add: No.1, Haobang Road, Yulin Economic Development Zone, Guangxi, China,
 Tel: +86-775-5802216 5802120 13877577356
 Fax: +86-775-5802213
 E-mail: HOTBORN@VIP.163.COM
 Web: WWW.HOTBORN.COM