

Professional Manufacture For Imaging Table.  
Best Quality In China.



# PRODUCT CATALOGUE



SUPPLY TABLES TO  
**GE, SIEMENS, PHILIPS, ZIEHM**  
AND KEEP A LONG TERM COOPERATION.



## BT-50E(II) SPINAL SURGERY TABLE

### POSITIONING CAPABILITIES

Prone  
Supine  
Lateral

Spinal Surgery Positioning with Exceptional C-arm, and compatible with Intelligent Surgical robot.

Follow the application of the Navigation system in Orthopedics operation, the Intelligent Orthopedics Surgery Robot is extensively used. Professional Spinal Surgery Table and precision positioning system is request for Intelligent center of orthopedics operation. This table is the good choice with precision positioning system and surgery navigation system.



## FEATURE

- ▶ Carbon fiber frame provides complete radiolucency.
- ▶ Compatible with variety of C-arm and O-arm
- ▶ Customized carbon fibre frame enhance positioning patients with varying body types.
- ▶ Cervical Traction Vector maintains the alignment of the traction rope.
- ▶ Multifunctional head rest is good perform on the flexion and head deflection of the cervical spine .
- ▶ The open base offers improved ergonomics and surgeon access.
- ▶ High-quality big casters and braking system allow clinicians to
- ▶ easily move the table.
- ▶ Micro-motors and accurate slide ensure quick, smooth, and balance table movement.
- ▶ 304 stainless steel design for durability and easy cleaning
- ▶ High-quality material, advanced and mature technology ensure the quality.





**COMPATIBLE WITH INTRA-OPERATIVE 3D IMAGING SYSTEM AND INTELLIGENT SURGICAL ROBOT.**

## COMMON PROCEDURES

- Spinal micro-operation positioning .
- Anterior/Posterior Fusions– cervical to sacrum
- Surgical Correction of Deformities
- Laminectomies
- Decompressions
- Kyphoplasties
- Osteotomies
- Spinal Fusions–Lateral, Anterior, Posterior, Transforaminal
- Intervertebral disc electrotherapy

## SPECIFICATIONS

- ▶ Table top(L\*W): 2000\*500 mm
- ▶ Height: 650~950 mm
- ▶ Trendelenburg & Reverse~: 10°
- ▶ Lateral tilt: 19°
- ▶ Load capacity: 200 kg
- ▶ Power: 220V 50Hz, 350VA



## STANDARD ACCESSORIES

- ▶ Mattress- 1 set
- ▶ Cable controller-1 set
- ▶ Buttack support frame (S)-1 pair
- ▶ Buttack support frame (M)-1 pair
- ▶ Buttack support frame (L)-1 pair
- ▶ Dual chest support frame (with pad)-1 set
- ▶ Knelt prone bag-1 unit
- ▶ Buttack strap-2 pieces
- ▶ Long rail clamp -2 sets
- ▶ Short rail clamp-2 sets
- ▶ Shoulder & Lateral support-1 pair
- ▶ Fixed strap-1 set
- ▶ IV Pole- 1 set
- ▶ Anaesthesia screen-1 set
- ▶ Head rest(FJ25)-1 set
- ▶ Arm board(FJ34)-1 pair
- ▶ Fixed strap (FJ07)-1 pair