

# DC-6000A

## CATHETERIZATION TABLE



TRANSLUCENT LENGTH : 1900mm

GOOD PERFORMANCE FOR VARIOUS C-ARM

## DATA



**HOTBORN**

Table top:	2450mm×490/635mm×40mm
Height (Electric):	755 ~ 1115mm
Vertical movement(Electric):	360mm
Length of table top without metal:	1178mm
Translucent Length:	1900mm
translucent width:	490mm
Head end dangling length:	2050mm
Foot end dangling length:	1225mm
Longitudinal sliding:	1200 mm
Lateral sliding:	310 mm
Load:	250kg



## PERFORMANCE



- Control system: adopt the precision Micro-motor, which is gear mesh system, to ensure the stable performance.
- Table top is made with good quality carbon fibre, whole piece without metal.
- Movement: vertical movement is by electro-mechanical drive, table top sliding is floating movement, is convenient and easy control.
- Center control foot brake
- Good quality wheel, integral forming, stable and wear proof.
- 304 Stainless steel base cover, low base to ensure good performance for C-arm.
- This table is special design for kinds of mobile Middle、small C-arm, for the intervention operations such as systemic vascular surgery, ERCP, Liver, Tumor, Vasculitis, pain etc. Product in good technology and quality.

## CONTENT

Main machine	1 set
Carbon fibre table top	1 set
Mattress	1 set
Hand controller	1 set
Arm support	1 set
Infusion lop	1 set

## HOTBORN

Add: No.8 Keyuan Dong Wu Road , Nanning New& Hi-Tech Industrial  
Development Zone, Nanning , Guangxi , China  
Tel: +86-771-3947312 13877577356  
Fax: +86-771-3947312  
E-mail: NNHOTBO@VIP.163.COM  
Web: WWW.HOTBORN.COM

Add: No.1, Haobang Road, Yulin Economic Development Zone, Guangxi, China,  
Tel: +86-775-5802216 5802120 13877577356  
Fax: +86-775-5802213  
E-mail: HOTBORN@VIP.163.COM  
Web: WWW.HOTBORN.COM

Note: The specifications are subject to change without prior notice for technical development